Julabo
THE TEMPERATURE CONTROL COMPANY

CORIO Series | DYNEO Series | HighTech Series

- Controle de temperatura preciso
- Mostradores luminosos brancos, bem visíveis à distância
- Operação simples



CORIO Series | **DYNEO Series** | HighTech Series

- Powerful circulation pumps – infinitely adjustable electronically
- Easy to switch between internal and external circulation
- Electronic fill level display
- Integrated programmer
- Multilingual display with intuitive user guidance
- Integrated Pt100 external sensor connection
- RS232 and anlog interfaces optional



CORIO Series | DYNEO Series | **HighTech Series**

- Elevadas potências de aquecimento
- Bombas de agitação potentes - ajustáveis eletronicamente em níveis
- Funções inteligentes para alarme e proteção – para maior segurança
- Sistema de aviso antecipado singular para nível baixo (DBGM)
- Interfaces digitais e analógicas – para uma comunicação flexível
- Monitoramento e operação sem fio através de WirelessTEMP (acessório)

